


DEQAS April 2017 - 1,25-dihydroxyvitamin D Method Means * (+/- 1 SD)



	Sample 401	Sample 402	Sample 403	Sample 404	Sample 405
■ DiaSorin Liaison XL (n=91)	70.9	113.2	135.2	93.2	162.3
■ DIASource CT (n=5)	71.2	121.8	151.2	100.8	151.3
■ IDS EIA (n=5)	85.4	154.1	182.1	120.3	201.9
■ IDS RIA (n=15)	100.5	175.3	199.7	134.3	202.0
■ IDS iSYS (n=18)	88.5	146.2	183.3	113.6	178.9
■ IDS iSYS New (n=12)	71.6	120.7	164.4	91.1	158.5
■ LC-MS/MS (n=13)	76.3	128.1	163.3	99.6	177.4
■ ALTM (n=160)	75.6	123.6	150.3	99.6	167.2