


## DEQAS April 2018 - 1,25-dihydroxyvitamin D Method Means \* ( $\pm$ 1 SD)



	Sample 421	Sample 422	Sample 423	Sample 424	Sample 425
■ DiaSorin Liaison XL (n=96)	62.2	79.1	153.9	132.0	100.7
■ DIASource CT (n=5)	77.5	83.1	175.9	152.6	112.9
■ IDS EIA (n=6)	61.5	73.1	148.1	123.5	101.4
■ IDS RIA (n=12)	74.1	82.7	194.1	157.4	123.6
■ IDS iSYS (n=10)	71.0	81.1	205.6	162.0	121.8
■ IDS iSYS New (n=20)	54.9	62.1	165.6	122.2	102.7
■ LC-MS/MS (n=12)	71.6	84.2	145.9	119.3	96.4
■ ALTM (n=161)	63.7	77.6	159.8	132.9	103.3

\* Methods with 5 or more results returned